

A close-up photograph of a person wearing a white lab coat and blue nitrile gloves. They are holding a red cylindrical pressure reactor with a stainless steel base and various fittings. A pressure gauge is attached to the top of the reactor, showing a reading of approximately 1000 psi. The background is a blurred laboratory setting.

# Training Pressure Reactors

Eric Hwang | YB Coporation. South Korea



BERGHOF

Products +  
Instruments


1.1

---

# Product Portfolio Overview

# Product Family

DAB



DB



BR | NR



Flexibility, Accessories, Functionality, Volume




Handling, Corrosion Resistance, robustness

# Range of Application

## Process sizes:

- Laboratory reactors (from 10 mL up to 2 L)
- Kilo Lab reactors (from 1 L up to 20 L)
- Pilot Plant reactors (from 5 L up to 200 L)
- Production reactors (from 50 L up to >1000 L)



# 1.2

---

## Product Portfolio

### Digestec DAB Vessels

# Simple Pressure Reactor?

# Digestec DAB-2/3

200 bar

250 °C

25 – 420 mL

3 vessel sizes

- High-quality stainless steel pressure vessels
- TFM-inserts, totally metal free
- Each pressure vessel is pre-examined by the TÜV
- Rupture discs as safety device
- Controlled ventilation of vessels during opening
- Simple handling
- Conventional heating via heating blocks
- Low costs



DAB-3



DAB-2

# Design of Digestec Pressure Vessels DAB-2/3



200 bar

250 °C

25 – 420 mL

3 vessel sizes

1. Pressure vessel
2. Bayonett quick closure
3. Base plate
4. Rupture disc holder
5. Disc springs
6. PTFE liner with cap




# Accessories of Digestec Pressure Vessels DAB-2/3

## Mounting tool for DAB-2/3 (p/n 5218060)


1. Assembly station
2. Torque wrench
3. Socket
4. PTFE paste
5. Brush
6. Mandrel for exchanging of the disc springs
7. Face spanner wrench



# DAH Heating Systems for Digestec Pressure Vessels DAB-2/3



1. Heating block
2. Power socket
3. Adjustable overtemperature protection
4. Thermocouple
5. Controller



DAH-902 with BTC-1000

# 1.3

---

## Product Portfolio

### Highpreactor DB Series

# But What if the Customer Wants to Measure $p$ and $T$ Inside the Reactor? → DB Series

200 bar

230 °C

300 – 2000 mL

6 vessel sizes

- High-quality stainless steel pressure vessels
- Complete PTFE lining
- Each reactor is build according to PED 2014/68/EU
- Rupture discs as safety device
- Pressure and temperature sensor available
- Stirring via magnetic stirring bar and hot plate (DB-300/500)
- Heating via double jacket (BTM) and/or electrical heating device (BHM) - see BR reactors



DB-500

DB-1000

# Constructing Materials / Sealing

## Liners made from PTFE for optimal chemical resistance


200 bar

230 °C

300 – 2000 mL

6 vessel sizes

- Removable PTFE inserts
  - Hermetically sealed between lid liner and insert
  - Dip tubes for thermocouple made from PTFE
  - PTFE sealing
  - $T$  up to 230 °C
  - $p$  up to 200 bar
- all parts in contact with gas and liquid phase are protected against chemical attack



# Design of Highpreactor DB Series

200 bar

230 °C


300 – 2000 mL

6 vessel sizes

1. Reactor vessel
2. Closing clamp
3. Opening tool
4. Reactor lid
5. Rupture disc holder
6. Vent valve
7. Dip tube for thermocouple
8. Pressure sensor

The following yellow parts inside the reactor are made of PTFE

1. Insert
2. O-ring sealing
3. Lid lining
4. Lining of the dip tube
5. Lining of the connections



# 1.4

---

## Product Portfolio

### Highpreactor BR Series

# BR Reactors with and w/o PTFE Insert

200 bar

230 / 300 °C

25 – 5700 mL

11 vessel sizes

- High-quality stainless steel pressure vessels
- Each reactor is build according to PED 2014/68/EU
- **PTFE insert and PTFE lid lining**
- Rupture discs as safety device
- O-ring as sealing system
- **Quick closing system – no wrench needed**
- Various armatures, configuration acc. customer request
- Magnetic coupled overhead stirring with different stirrer geometries
- Heating via double jacket (BTM) and/or electrical heating device (BHM)



BR-500



BR-2000

# Heating devices

200 bar

230 / 300 °C

25 – 5700 mL

11 vessel sizes

## Two alternative heating devices

### 1. BHM „Berghof Heating Mantle“

- Electrical heating device with over temperature protection
- Controller BTC/BDL needed

### 2. BTM „Berghof Thermostated Mantle“

- Double jacket for heating with oil / heating fluid
- Heating circulator needed



BHM



BTM

# Overview BR Reactor Series

**BR-25/40**



**BR-100/200**



**BR-300/500/700**



**BR-800/1000/1500/2000**



**BR-4000**



Accessories, Functionality, Volume, Diameter of Vessel and Lid

# General Design of Highpreactor BR Series

200 bar

230 °C


25 – 5700 mL

11 vessel sizes

1. Reactor vessel
2. Closing clamp
3. Opening tool
4. Reactor lid
5. Rupture disc holder
6. Vent valve
7. Dip tube for thermocouple
8. Manometer / pressure gauge
9. Magnetic clutch

The following yellow parts inside the reactor are made of PTFE if needed

1. Insert
2. O-ring sealing
3. Lid lining
4. Lining of the dip tube
5. Sampling tube
6. Stirrer



# Constructing Materials / Sealing

200 bar

230 / 300 °C

25 – 5700 mL


11 vessel sizes

## Two alternative sealing positions at BR-300 and bigger

1. Outer O-ring for working with stainless steel reactors only
2. Inner O-ring for working with PTFE-lined reactors

## Different materials of construction


- Stainless steel 316Ti (1.4571) as standard material
- Hastelloy C-22 (2.4602)
- Other materials and coatings, e.g. Tantalum, on request



# BR-25/40 Reactors with and w/o PTFE Insert

	200 bar	230 / 300 °C	25 – 45 mL
Series	Volume vessel (metal-version)		Volume PTFE insert
BR-25	30 mL		25 mL
BR-40	45 mL		40 mL


- Attention! Different vessels for reactors with or w/o PTFE
- Attention! With PTFE insert → 230 °C max
- Reactor lid as screw cap
- Stirring via magnetic stirring bar, no magnetic clutch available
- Heating via laboratory heating plate



# BR-100/200 Reactors with and w/o PTFE Insert

	200 bar	230 / 300 °C	25 – 200 mL
Series	Volume vessel (metal-version)	Volume PTFE insert	
BR-100	100 mL	50/75/100/150 mL	
BR-200	200 mL	200 mL	

- Attention! Different vessels for reactors with or w/o PTFE
- Attention! With PTFE insert → 230 °C max
- Closing clamp
- Stirring via magnetic stirring bar, no magnetic clutch available
- Heating via laboratory heating plate



# BR-300...2000 Reactors with and w/o PTFE Insert


200 bar

230 / 300 °C

300 – 2200 mL

Series	Volume vessel	Volume PTFE insert
BR-300	400 mL	300 mL
BR-500	600 mL	500 mL
BR-700	900 mL	700 mL
BR-1000	1200 mL	1000 mL
BR-1500	1700 mL	1500 mL
BR-2000	2200 mL	2000 mL

- Magnetic coupled overhead stirring with bottom bearing
- Electric heating mantle (BHM) with controller and/or thermostated mantle (BTM)
- **Attention! BR-300 without bottom bearing**



# BR-4000 Reactors with and w/o PTFE Insert

150 bar

230 / 300 °C

4000 / 5700 mL

Series

Volume vessel

Volume PTFE insert

BR-4000

5700 mL

4000 mL

- Magnetic coupled overhead stirring with bottom bearing
- Electric heating mantle (BHM) with controller and/or thermostated mantle (BTM)
- Closing clamp
- Lift system or laboratory cart for easy handling



# 1.5

---

## Product Portfolio

### Highpreactor NR Series

# NR Reactors with and w/o PTFE Insert



25 bar

230 / 300 °C

1500 – 11400 mL

Serie	Volume vessel	Volume PTFE insert
NR-1500	1800 mL	1500 mL
NR-3000	3500 mL	3000 mL
NR-6000	6900 mL	6000 mL
NR-10000	11400 mL	-

- Available with and w/o PTFE insert, except NR-10000
- Bottom drain valve is also available, but only w/o PTFE
- Magnetic coupled overhead stirring with bottom bearing
- Electrical heating mantle with controller or isolated double jacket, directly welded onto the reactor vessel
- Customization in size, materials and features possible



# 2.1

---

## Applications

DAB, DB, NR and BR Systems

# Examples of Applications

## DAB Vessels


**Material testing / Aging tests / Exposure tests**

Specimens are treated in liquids / gases at elevated temperatures for defined time periods

Typical parameters:

- Specimen: Elastomers, Thermoplasts, Metal parts
- Liquids: cooling liquids, mineral oils, salt water, fuels
- $T = 150\text{ °C}$  and  $t$  up to 1000 h

Customers: **Companies like Daimler / Mercedes Benz, Continental, Trelleborg, Baker & Hughes, ExxonMobil...**



**Tensile test**



**DAB-3 XXL with specimen holder**

# Examples of Applications

## DAB Vessels

### Solvothermal synthesis

Synthesis of porous and non porous solids under solvothermal conditions / crystallization step

Typical parameters:

- Products: Zeolites, mesoporous metal oxides, nanoparticles...
- Solvents: Water, Acids, Bases, Ethanol...
- $T = 110$  to  $250$  °C and  $t$  up to 168 h

Customers: **Academic research, Wacker, BASF, Evonik, Bayer, ExxonMobil, Shell...**



DAB-3



DAB-2

# Examples of Applications

## DB Vessels

### Corrosion tests

Specimens are treated in liquids / gases at elevated temperatures for defined time periods

Typical parameters:

- Specimen: Metal coupons, metal parts...
- Liquids: cooling liquids, salt water, process water
- $T = 150\text{ °C}$  and  $t$  up to 672 h

Customers: **Companies like Linde, Kurita, Daimler / Mercedes Benz...**



DB-500

# Examples of Applications

## DB Vessels

### Digestion / Acid treatment

Digestion and/or acid treatment of rocks, metals, metal oxides, silicates...

Typical parameters:

- Liquids: concentrated inorganic acids / mixtures of acids
- $T$  up to 230 °C and  $t$  up to 48 h

Customers: **Academic research, Rusal, Daimler / Mercedes Benz, Fraunhofer Institute...**



DB-1000

# Examples of Applications

## BR-25/40 and BR-100/200

### Small scale synthesis

Organic and Inorganic small scale synthesis

Typical parameters:

- Reactions: Hydrogenations, Metathesis, Gas-liquid reactions...
- Solvents: Aqueous solutions, Alcohols, Alkanes, Aromatics...
- $T$  up to 280 °C,  $p$  up to 180 bar and  $t$  up to 6 h

Customers: **Industrial and Academic research, Pharma research, Boehringer Ingelheim, Roche Diagnostics...**



BR-25

BR-100

# Examples of Applications

## BR-300...2000

Organic and Inorganic Synthesis

Organic and Inorganic synthesis

Typical parameters:

- Reactions: Hydrogenations, Polymerizations, Gas-liquid reactions...
- Solvents: Aqueous solutions, Alcohols, Alkanes, Aromatics...
- $T$  up to 280 °C,  $p$  up to 150 bar and  $t$  up to 6 h

Customers: **Industrial and Academic research, Pharma research, Continental, Panacol, Wacker, Dow, Evonik, Bayer, BASF, Clariant...**



BR-500

# Examples of Applications

## BR-700 and 2000

**Material testing / Aging tests / Exposure tests**

Specimens are treated in liquids / gases at elevated temperatures for defined time periods

Typical parameters:

- Specimen: Elastomers, Thermoplasts, Metal parts
- Liquids: cooling liquids, synthetic oils, fuels
- $T = 150\text{ °C}$  and  $t$  up to 1000 h

Customers: **Companies like Daimler / Mercedes Benz, Continental, Trelleborg, Baker & Hughes, ExxonMobil, Freudenberg, Vaxxilon, Schaeffler...**



Sample holder for BR-700 and BR-2000

# Examples of Applications

## BR-4000

### Waste treatment

Sub-process of the treatment of industrial waste for recycling of precious metals

Typical parameters:

- Educts: cyanide solutions with precious metal ions
- Liquids: concentrated alkaline solutions like caustic soda and potash
- $T = 200\text{ }^{\circ}\text{C}$  and  $t$  up to 24 h

Customers: **Recycling companies like Remondis, Alba...**



Electronic waste material  
before crushing / milling and  
treatment with cyanide solution

# Examples of Applications

## BR-4000

### Hydrothermal Carbonization

Conversion of biomass to brown coal

Typical parameters:

- Feedstock: All kinds of biomass
- Liquids: slightly acid solutions
- $T = 230\text{ }^{\circ}\text{C}$  and  $t$  up to 72 h at hydrothermal conditions

Customers: **Academic research and companies like Spina Group, Suncoal Industries...**



BR-4000

Thank you | 감사합니다



Your Best Laboratory

[berghof-instruments.com](http://berghof-instruments.com)



Products +  
Instruments