

SHP Steriltechnik AG



Schloss Detzel 1
39345 Detzel Schloss / Satuelle
Tel.: 039058 / 9762-0
Fax: 039058 / 9762-22



Application....Solution

Sterile – what does that mean?

- ◆ A subject or a liquid is sterile if it is free of augmentable microorganism.
- ◆ sterile does NOT mean clean
- ◆ that is to say a sterilization process is no cleaning
- ◆ A dirty glas remains dirty after sterilization.
 - it is sterile
 - thus free of augmentable microorganism
 - But not clean!

Sterile – what does that mean?

Air – the disturbing factor of sterilizing

Where is air, there is **no** steam!

Where is no steam, there is **no** adequate temperature!

Where is **no** adequate temperature, there is **no successful sterilization** possible.




Sterile – what does that mean?

An efficient evacuation of the air out of the sterilizing chamber is the basis for a successful sterilization!


Methods:

- ◆ gravitation
- ◆ pressure purge
- ◆ de-aeration by pre-vacuum
 - simple pre-vacuum
 - fractionated pre-vacuum


Sterilizing goods



Sterilizing goods



Sterilizing goods




Sterilizing goods

mixed fractions

mixture of:
liquids and solids
(waste sterilization)


Sterilizing goods




Methods of de-aeration for liquids

Gravitation


- ◆ chamber: filled with **air** and loading
- ◆ **steam** reaches the chamber, pressure increases
- ◆ presses **air** down
- ◆ outlet valve is opened
- ◆ **air** escapes
- ◆ at 96°C post de-aeration for 10min




Re-cooling systems for liquids



air re-cooling




water re-cooling by
cooling coil



water re-cooling by
double jacket

Re-cooling systems for liquids


Air re-cooling



- ◆ ambient air is channeled on the surface of the chamber
- ◆ cooling only where air hits the chamber
 - ineffective
 - slow
 - small area of contact
- ◆ no cooling with counter pressure possible

Re-cooling systems for liquids


Water re-cooling by cooling coil



- ◆ demineralized water goes through the cooling coil
- ◆ cooling only where cooling coil touches the chamber
 - small contact area
- ◆ cooling with counter pressure possible

Re-cooling systems for liquids

Water re-cooling by double wall



- ◆ demineralized water goes through the complete double jacket
- ◆ double jacket encloses the chamber
 - big area of contact
 - effective and fast cooling
- ◆ cooling with counter pressure possible

Re-cooling systems for liquids

Water re-cooling by double wall

The bigger the area of contact

- the bigger is heat exchange
- the faster and efficient is re-cooling

Conclusion:

**double wall = most efficient
re-cooling system**



Sterilizing goods

Liquids – open or lightly closed vessels


- ◆ without active re-cooling
 - ◆ programs:
 - Liquids SL (self re-cooling)
 - Liquids ST (clocked steam outlet)
- Laboklav 25B
- Laboklav 55-195B



Sterilising goods

Liquids – open or lightly closed vessels

Liquid ST
(without active re-cooling)



- ◆ de-aeration by gravitation
- ◆ slow and passive re-cooling (no re-cooling system)
 - only slow pressure purge

Laboklav 25B

Laboklav 55-195B

Sterilizing goods

Liquids – open or lightly closed vessels

- ◆ active re-cooling
 - loss of liquids possible
 - ◆ program
 - Liquids RM
- Laboklav 25M
 - air re-cooling
- Laboklav 55-195M
 - water re-cooling by double jacket




Sterilizing goods

Liquids – open or lightly closed vessels

Liquids RM

(Liquid programm with re-cooling by double jacket)



- ◆ de-aeration by gravitation
- ◆ re-cooling by double jacket

- loss of liquids in open vessels app. 6 – 12%
- boiling retardation possible
- concentration in liquid changes

Laboklav 25M (air re-cooling)

Laboklav 55-195M

Sterilizing goods

Liquids – open or lightly closed vessels

- ◆ active re-cooling and counter pressure
- ◆ but slower than just re-cooling by double jacket
- ◆ program
 - Liquids RO
- ◆ Laboklav 25
 - not possible
- Laboklav 55-195MS




Sterilizing goods

Liquids – open or lightly closed vessels

Liquids RO

(with re-cooling by double jacket & counter



- ◆ de-aeration by gravitation
- ◆ re-cooling by jacket with counter pressure

- no loss of liquids
- no boiling retardation
- but slower than just re-cooling by jacket
- temperature gradients possible

Laboklav 55-195MS

Sterilizing goods

Liquids – open or lightly closed vessels

- active re-cooling, counter pressure and fan inside chamber
 - no temperature gradients inside chamber
 - program
 - Liquids RO
- Laboklav 25
→ not possible
→ Laboklav 55-195MSL




Sterilizing goods

Liquids – open or lightly closed vessels

Liquids RO

(with re-cooling by double jacket, counter



◆ de-aeration by gravitation

◆ re-cooling by jacket with counter pressure

- no loss of liquids
- no boiling retardation
- as fast as re-cooling by jacket
- no temperature gradients

Laboklav 55-195MSL

Sterilizing goods

Liquids – pressure tight closed vessels

- ◆ water re-cooling, counter pressure and fan inside chamber
- ◆ program
 - Liquids RG
- ◆ Laboklav 25
 - not possible
 - Laboklav 55-195MSL




Sterilizing goods

Liquids – pressure tight closed vessels

Liquid RG

(with re-cooling by jacket & counter pressure, already during preheating, fan)



Laboklav 55-195MSL

- ◆ no real de-aeration
 - vibrant de-aeration
- ◆ heating with counter pressure
 - overpressure in the vessels because of fast increasing temperature)
 - counter pressure = backpressure
- ◆ counter pressure during sterilization phase
- ◆ jacket re-cooling with counter pressure and fan inside chamber
 - no temperature gradients inside chamber
 - shorter cooling time


Please consider on units with counter pressure:

Laboklav 55-195 MS, MSL, MSV & MSLV

- ◆ autoclaves with the option **counter pressure** do not work properly without counter pressure
- ◆ compressed air is needed for each program
 - compressed air is used for aeration of the chamber
 - air which is led into chamber is pressed through the sterile filter by the compressor and goes afterwards into the chamber

**without aeration → no equalisation
→ unit cannot be opened!**


Sterilizing goods



Methods of de-aeration for solids

Pressure purge

- ◆ for easy to de-aerate solids
- ◆ **steam** injecting into chamber
 - pressure rises, outlet valve closed
- ◆ outlet valve opens, **air** removes
- ◆ meanwhile no injecting of **steam**
- ◆ process is repeated up to 6x



Sterilizing goods

Solids – glasware and instruments


- ◆ de-aeration by pressure purge
 - ◆ without drying in vacuum
 - ◆ program
 - Instruments ST
- Laboklav 25B
- Laboklav 55-195B



Sterilizing goods

Solids – glasware and instruments

Instruments ST (solid programm)



- ◆ de-aeration by pressure purge
- ◆ no drying in vacuum


Laboklav 25B

Laboklav 55-195B

Methods of de-aeration for solids

Fractionated pre-vacuum

- ◆ method of de-aeration for difficult to de-aerate solids
- ◆ simple pre-vacuum
 - until adjusted setpoint
- ◆ or fractionated pre-vacuum
 - 3 fractions of de-aeration
- ◆ change of steam injecting and vacuum



Sterilizing goods

Solids – hollow items and hoses

- ◆ vacuum required
 - simple or fractionated pre-vacuum
- ◆ drying vacuum recommendable
- ◆ program
 - Instruments FV
 - Instruments FV with VT



→ Laboklav 25V


→ Laboklav 55-195V

Sterilizing goods

Solids – hoses and hollow items

Instruments FV

(Solid programm with fractionated pre-vacuum)



- ◆ de-aeration by fractionated pre-vacuum
- ◆ no drying vacuum
 - e.g. for pipette tips
 - but necessary for hoses (must be sterile inside)


Laboklav 25V,55-195V

Sterilizing goods

Solids – hoses and hollow items

Instruments FV with VT

(Solid programm with pre-vacuum & drying vacuum)



- ◆ de-aeration by fractionated pre-vacuum
- ◆ drying vacuum

- heating of the chamber
- condensate on the sterilising good becomes steam
- vacuum pump pulls steam out of chamber

Laboklav 25V,55-195V

Sterilizing goods

Solids – textiles and porous goods

- ◆ de-aeration by vacuum
 - fractionated pre-vacuum
- ◆ drying in vacuum required
- ◆ program
 - Instruments FV with VT

→ Laboklav 25V

→ Laboklav 55-195V




Sterilizing goods

Solids – hoses and hollow items

Instruments FV with VT

(Solid programm with pre-vacuum & drying vacuum)




- ◆ de-aeration by fractionated pre-vacuum
- ◆ drying vacuum

- heating of the chamber
- condensate on the sterilising good becomes steam
- vacuum pump pulls steam out of chamber

Laboklav 25V,55-195V

Sterilization goods

Waste strilization



Sterilizing goods

Mixed fractions – solids and amounts of liquids <0,5l

- ◆ in open waste bags
 - no pre-vacuum required
 - fill in 0,5l demineralized water into opened waste bag
- ◆ program
 - Destruct ST

→ Laboklav 25B

→ Laboklav 55-195B




Sterilizing goods

Mixed fractions – solids and amounts of liquids <0,5l

Destruct ST
(solid program)

- ◆ de-aeration by pressure purge
- ◆ no drying in vacuum



Laboklav 25B

Laboklav 55-195B

Sterilizing goods

Mixed fractions – solids and amounts of liquids <0,5l

- ◆ in closed waste bags
 - fractionated pre-vacuum required
- ◆ program
 - Destruct FV

→ Laboklav 25V

→ Laboklav 55-195V




Sterilizing goods

Mixed fractions – solids and amounts of liquids <0,5l

Destruct FV

(solid program with fractionated pre-vacuum)



- ◆ de-aeration by fractionated pre-vacuum
- ◆ no drying in vacuum

Laboklav 25V

Laboklav 55-195V

Sterilizing goods

Mixed fractions – solids and amounts of liquids $>0,5l$

- ◆ basic unit is adequate
 - ◆ the bigger the amount of liquids, the more recommendable is a re-cooling system
 - ◆ program
 - Liquid waste
 - or all other liquid programs
- Laboklav 25B & M
- Laboklav 55-195B, M, MS, MSL




Sterilizing goods

Mixed fractions – solids and amounts of liquids >0,5l

Liquid waste
(liquid program)

- ◆ de-aeration by gravitation
- ◆ passive re-cooling



- re-cooling system shorten batch time
- the more liquids are sterilised, the more recommendable is a re-cooling system

Laboklav 25B, M

Laboklav 55-195B, M, MS, MSL

Sterilization goods

Big amount of liquid

capacity up to
100 litres of liquid



Laboklav Wastewater

Thank you for your attention!